



Material			
No.	Plan	Article	Quantity
1	X003872	1901134	1
2	Y008382	1687104 (L)	1
	Y013417	1687204 (R)	1
3	Y012555	1676204	1
4	Y012556	1676100	1
5	---	1418700	1
6	M10 DIN934	1001672	1
7*	M6 DIN985	1001680	3
8	M8 DIN125	1003199	1
9	M6x20 ISO7380	1006999	1
11*	M6x30 ISO7380	1009588	2
12	M10x19	1991432	1
13**	X003698	1004550	2
14**	M5 DIN985	2994006	2
15	M6x25 ISO7380	1006998	1
16**	M5x20 ISO7380	1998368	2

* For Eclips / 708D

** For BI5 V500 / Trigo / Turios

Supported models:
Eclips / 708D / V500 / Trigo / Turios

19033364 / 19033365
1913829 / 1913830

E000523 / H

B66-O ADJUSTABLE ARMREST B66-O VERSTELBAAR				Plan No. E000523 / H
Scale: 1:2	Drafted	Verified	Confirmed	
Date	28/10/2008	09/11/2022	09/11/2022	
Initials	SD	SD	PJ	